



TRPV2 Antibody, HRP conjugated

Product Code	CSB-PA897571LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5S1
Immunogen	Recombinant Human Transient receptor potential cation channel subfamily V member 2 protein (1-117AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transient receptor potential cation channel subfamily V member 2 (TrpV2) (Osm-9-like TRP channel 2) (OTRPC2) (Vanilloid receptor-like protein 1) (VRL-1), TRPV2, VRL
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TRPV2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.