




ABL2 Antibody

Product Code	CSB-PA001106LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42684
Immunogen	Recombinant Human Tyrosine-protein kinase ABL2 protein (862-1064AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Abelson tyrosine-protein kinase 2 (EC 2.7.10.2) (Abelson murine leukemia viral oncogene homolog 2) (Abelson-related gene protein) (Tyrosine-protein kinase ARG), ABL2, ABL2 ARG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ABL2

Image



Western Blot

Positive WB detected in: U87 whole cell lysate, U251 whole cell lysate

All lanes: ABL2 antibody at 4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 129, 125, 127, 115, 118,

116, 117, 61 kDa

Observed band size: 115 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.