



ABL2 Antibody, Biotin conjugated

Product Code	CSB-PA001106LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42684
Immunogen	Recombinant Human Tyrosine-protein kinase ABL2 protein (862-1064AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Abelson tyrosine-protein kinase 2 (EC 2.7.10.2) (Abelson murine leukemia viral oncogene homolog 2) (Abelson-related gene protein) (Tyrosine-protein kinase ARG), ABL2, ABL2 ARG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ABL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.