



# ABL2 Antibody, Biotin conjugated

Product Code	CSB-PA001106LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42684
Immunogen	Recombinant Human Tyrosine-protein kinase ABL2 protein (862-1064AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Abelson tyrosine-protein kinase 2 (EC 2.7.10.2) (Abelson murine leukemia viral oncogene homolog 2) (Abelson-related gene protein) (Tyrosine-protein kinase ARG), ABL2, ABLL ARG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ABL2