



ADAMTS13 Antibody, HRP conjugated

Product Code	CSB-PA001301LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q76LX8
Immunogen	Recombinant Human A disintegrin and metalloproteinase with thrombospondin motifs 13 protein (738-905AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	A disintegrin and metalloproteinase with thrombospondin motifs 13 (ADAM-TS 13) (ADAM-TS13) (ADAMTS-13) (EC 3.4.24.87) (von Willebrand factor-cleaving protease) (vWF-CP) (vWF-cleaving protease), ADAMTS13, C9orf8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ADAMTS13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.