



AGPS Antibody, Biotin conjugated

Product Code	CSB-PA001457LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00116
Immunogen	Recombinant Human Alkyldihydroxyacetonephosphate synthase, peroxisomal protein (279-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alkyldihydroxyacetonephosphate synthase, peroxisomal (Alkyl-DHAP synthase) (EC 2.5.1.26) (Aging-associated gene 5 protein) (Alkylglycerone-phosphate synthase), AGPS
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	AGPS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.