



C4A Antibody, Biotin conjugated

Product Code	CSB-PA003949LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C0L4
Immunogen	Recombinant Human Complement C4-A protein (1027-1186AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complement C4-A (Acidic complement C4) (C3 and PZP-like alpha-2-macroglobulin domain-containing protein 2) [Cleaved into: Complement C4 beta chain; Complement C4-A alpha chain; C4a anaphylatoxin; C4b-A; C4d-A; Complement C4 gamma chain], C4A, CO4 CPAMD2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	C4A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.