



CALR Antibody, Biotin conjugated

Product Code	CSB-PA004458YD01PI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28491
Immunogen	Recombinant Pig Calreticulin protein (18-417AA)
Raised In	Rabbit
Species Reactivity	Pig
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calreticulin (CRP55) (Calregulin) (Endoplasmic reticulum resident protein 60) (ERp60) (HACBP), CALR
Immunogen Species	Sus scrofa (Pig)
Research Area	Tags & Cell Markers
Target Names	CALR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.