



CDR2L Antibody, Biotin conjugated

Product Code	CSB-PA005104LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86X02
Immunogen	Recombinant Human Cerebellar degeneration-related protein 2-like protein (376-465AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cerebellar degeneration-related protein 2-like (Paraneoplastic 62 kDa antigen), CDR2L, HUMPPA
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CDR2L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.