



CEACAM16 Antibody, Biotin conjugated

Product Code	CSB-PA005158LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2WEN9
Immunogen	Recombinant Human Carcinoembryonic antigen-related cell adhesion molecule 16 protein (323-414AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carcinoembryonic antigen-related cell adhesion molecule 16 (Carcinoembryonic antigen-like 2), CEACAM16, CEAL2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CEACAM16
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.