



CLCNKB Antibody, Biotin conjugated

Product Code	CSB-PA005488LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51801
Immunogen	Recombinant Human Chloride channel protein CIC-Kb protein (538-687AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Chloride channel protein CIC-Kb (Chloride channel Kb) (CIC-K2), CLCNKB
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CLCNKB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.