




CPE Antibody

Product Code	CSB-PA005886LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16870
Immunogen	Recombinant Human Carboxypeptidase E protein (252-363AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carboxypeptidase E (CPE) (EC 3.4.17.10) (Carboxypeptidase H) (CPH) (Enkephalin convertase) (Prohormone-processing carboxypeptidase), CPE
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CPE

Image



Western Blot

Positive WB detected in: A375 whole cell lysate

All lanes: CPE antibody at 2.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 54, 50 kDa

Observed band size: 54 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.