



CPS1 Antibody, HRP conjugated

Product Code	CSB-PA005913LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31327
Immunogen	Recombinant Human Carbamoyl-phosphate synthase [ammonia], mitochondrial protein (1361-1500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carbamoyl-phosphate synthase [ammonia], mitochondrial (EC 6.3.4.16) (Carbamoyl-phosphate synthetase I) (CPSase I), CPS1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	CPS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.