



CRLF1 Antibody, HRP conjugated

Product Code	CSB-PA005982LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75462
Immunogen	Recombinant Human Cytokine receptor-like factor 1 protein (318-422AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytokine receptor-like factor 1 (Cytokine-like factor 1) (CLF-1) (ZcytoR5), CRLF1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CRLF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.