



DPPA5 Antibody, FITC conjugated

Product Code	CSB-PA007147LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6NC42
Immunogen	Recombinant Human Developmental pluripotency-associated 5 protein (1-116AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Developmental pluripotency-associated 5 protein (hDPPA5) (Embryonal stem cell-specific gene 1 protein) (ESG-1), DPPA5, ESG1
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	DPPA5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.