



DSCAM Antibody, HRP conjugated

Product Code	CSB-PA007192LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60469
Immunogen	Recombinant Human Down syndrome cell adhesion molecule protein (1711-1867AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Down syndrome cell adhesion molecule (CHD2), DSCAM
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	DSCAM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.