



DSCAM Antibody, Biotin conjugated

Product Code	CSB-PA007192LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60469
Immunogen	Recombinant Human Down syndrome cell adhesion molecule protein (1711-1867AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Down syndrome cell adhesion molecule (CHD2), DSCAM
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	DSCAM