



ECM1 Antibody, FITC conjugated

Product Code	CSB-PA007383LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16610
Immunogen	Recombinant Human Extracellular matrix protein 1 protein (386-540AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Extracellular matrix protein 1 (Secretory component p85), ECM1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	ECM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.