



ESRRB Antibody, HRP conjugated

Product Code	CSB-PA007836LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95718
Immunogen	Recombinant Human Steroid hormone receptor ERR2 protein (432-502AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Steroid hormone receptor ERR2 (ERR beta-2) (Estrogen receptor-like 2) (Estrogen-related receptor beta) (ERR-beta) (Nuclear receptor subfamily 3 group B member 2), ESRRB, ERB2 ESRL2 NR3B2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ESRRB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.