



# ESRRB Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA007836LD01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O95718
<b>Immunogen</b>	Recombinant Human Steroid hormone receptor ERR2 protein (432-502AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Steroid hormone receptor ERR2 (ERR beta-2) (Estrogen receptor-like 2) (Estrogen-related receptor beta) (ERR-beta) (Nuclear receptor subfamily 3 group B member 2), ESRRB, ERB2 ESRL2 NR3B2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	ESRRB