



EYA4 Antibody, Biotin conjugated

Product Code	CSB-PA007909LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95677
Immunogen	Recombinant Human Eyes absent homolog 4 protein (153-224AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Eyes absent homolog 4 (EC 3.1.3.48), EYA4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	EYA4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.