



FCHO1 Antibody, HRP conjugated

Product Code	CSB-PA008546LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14526
Immunogen	Recombinant Human F-BAR domain only protein 1 protein (45-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-BAR domain only protein 1, FCHO1, KIAA0290
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FCHO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.