




FOXL2 Antibody

Product Code	CSB-PA008827LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P58012
Immunogen	Recombinant Human Forkhead box protein L2 protein (116-216AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Forkhead box protein L2, FOXL2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FOXL2

Image



Western Blot

Positive WB detected in: A549 whole cell lysate
 All lanes: FOXL2 antibody at 3.7µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 39 kDa
 Observed band size: 49 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.