



GBX2 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA009305LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P52951 |
| Immunogen | Recombinant Human Homeobox protein GBX-2 protein (65-248AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Homeobox protein GBX-2 (Gastrulation and brain-specific homeobox protein 2), GBX2 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | GBX2 |