



GBX2 Antibody, Biotin conjugated

Product Code	CSB-PA009305LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52951
Immunogen	Recombinant Human Homeobox protein GBX-2 protein (65-248AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein GBX-2 (Gastrulation and brain-specific homeobox protein 2), GBX2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GBX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.