



GEMIN4 Antibody, HRP conjugated

Product Code	CSB-PA009363LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57678
Immunogen	Recombinant Human Gem-associated protein 4 protein (741-848AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gem-associated protein 4 (Gemin-4) (Component of gems 4) (p97), GEMIN4
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	GEMIN4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.