



GMDS Antibody, FITC conjugated

Product Code	CSB-PA009569LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60547
Immunogen	Recombinant Human GDP-mannose 4,6 dehydratase protein (243-346AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GDP-mannose 4,6 dehydratase (EC 4.2.1.47) (GDP-D-mannose dehydratase) (GMD), GMDS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GMDS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.