



GREM1 Antibody, Biotin conjugated

Product Code	CSB-PA009892LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60565
Immunogen	Recombinant Human Gremlin-1 protein (18-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gremlin-1 (Cell proliferation-inducing gene 2 protein) (Cysteine knot superfamily 1, BMP antagonist 1) (DAN domain family member 2) (Down-regulated in Mos-transformed cells protein) (Increased in high glucose protein 2) (IHG-2), GREM1, CKTSF1B1 DAND2 DRM
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GREM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.