



GSS Antibody, FITC conjugated

Product Code	CSB-PA009969LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48637
Immunogen	Recombinant Human Glutathione synthetase protein (290-418AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutathione synthetase (GSH synthetase) (GSH-S) (EC 6.3.2.3) (Glutathione synthase), GSS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GSS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.