



GTPBP1 Antibody, Biotin conjugated

Product Code	CSB-PA010031LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00178
Immunogen	Recombinant Human GTP-binding protein 1 protein (562-669AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GTP-binding protein 1 (G-protein 1) (GP-1) (GP1), GTPBP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GTPBP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.