



HNRNPAB Antibody, HRP conjugated

Product Code	CSB-PA010604LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99729
Immunogen	Recombinant Human Heterogeneous nuclear ribonucleoprotein A/B protein (1-68AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heterogeneous nuclear ribonucleoprotein A/B (hnRNP A/B) (APOBEC1-binding protein 1) (ABBP-1), HNRNPAB, ABBP1 HNRNPAB
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HNRNPAB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.