



HPN Antibody, Biotin conjugated

Product Code	CSB-PA010704LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05981
Immunogen	Recombinant Human Serine protease hepsin protein (294-413AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine protease hepsin (EC 3.4.21.106) (Transmembrane protease serine 1) [Cleaved into: Serine protease hepsin non-catalytic chain; Serine protease hepsin catalytic chain], HPN, TMPRSS1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	HPN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.