



IL1RAPL1 Antibody, FITC conjugated

Product Code	CSB-PA011624LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZN1
Immunogen	Recombinant Human Interleukin-1 receptor accessory protein-like 1 protein (564-679AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-1 receptor accessory protein-like 1 (IL-1-RAPL-1) (IL-1RAPL-1) (IL1RAPL-1) (Oligophrenin-4) (Three immunoglobulin domain-containing IL-1 receptor-related 2) (TIGIRR-2) (X-linked interleukin-1 receptor accessory protein-like 1), IL1RAPL1, OPHN4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	IL1RAPL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.