



INPP5B Antibody, HRP conjugated

Product Code	CSB-PA011734LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32019
Immunogen	Recombinant Human Type II inositol 1,4,5-trisphosphate 5-phosphatase protein (78-178AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Type II inositol 1,4,5-trisphosphate 5-phosphatase (EC 3.1.3.36) (75 kDa inositol polyphosphate-5-phosphatase) (Phosphoinositide 5-phosphatase) (5PTase), INPP5B, OCRL2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	INPP5B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.