



KCNJ5 Antibody, Biotin conjugated

Product Code	CSB-PA012058LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48544
Immunogen	Recombinant Human G protein-activated inward rectifier potassium channel 4 protein (348-419AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	G protein-activated inward rectifier potassium channel 4 (GIRK-4) (Cardiac inward rectifier) (CIR) (Heart KATP channel) (Inward rectifier K(+) channel Kir3.4) (IRK-4) (KATP-1) (Potassium channel, inwardly rectifying subfamily J member 5), KCNJ5, GIRK4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNJ5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.