



IST1 Antibody, Biotin conjugated

Product Code	CSB-PA012171LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53990
Immunogen	Recombinant Human IST1 homolog protein (172-284AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	IST1 homolog (hIST1) (Putative MAPK-activating protein PM28), IST1, KIAA0174
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	IST1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.