



# LGMN Antibody, Biotin conjugated

Product Code	CSB-PA012903LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99538
Immunogen	Recombinant Human Legumain protein (199-312AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Legumain (EC 3.4.22.34) (Asparaginyl endopeptidase) (Protease, cysteine 1), LGMN, PRSC1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LGMN