



LRPPRC Antibody, Biotin conjugated

Product Code	CSB-PA013105LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42704
Immunogen	Recombinant Human Leucine-rich PPR motif-containing protein, mitochondrial protein (901-1036AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich PPR motif-containing protein, mitochondrial (130 kDa leucine-rich protein) (LRP 130) (GP130), LRPPRC, LRP130
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	LRPPRC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.