



LSAMP Antibody, HRP conjugated

Product Code	CSB-PA013197LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13449
Immunogen	Recombinant Human Limbic system-associated membrane protein (168-287AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Limbic system-associated membrane protein (LSAMP) (IgLON family member 3), LSAMP, IGLON3 LAMP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	LSAMP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.