



LSAMP Antibody, HRP conjugated

Product Code	CSB-PA013197LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13449
Immunogen	Recombinant Human Limbic system-associated membrane protein (168-287AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Limbic system-associated membrane protein (LSAMP) (IgLON family member 3), LSAMP, IGLON3 LAMP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	LSAMP