



MVD Antibody, Biotin conjugated

Product Code	CSB-PA015246LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53602
Immunogen	Recombinant Human Diphosphomevalonate decarboxylase protein (69-221AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Diphosphomevalonate decarboxylase (EC 4.1.1.33) (Mevalonate (diphospho)decarboxylase) (MDDase) (Mevalonate pyrophosphate decarboxylase), MVD, MPD
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	MVD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.