



MX2 Antibody, Biotin conjugated

Product Code	CSB-PA015250LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20592
Immunogen	Recombinant Human Interferon-induced GTP-binding protein Mx2 protein (548-688AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon-induced GTP-binding protein Mx2 (Interferon-regulated resistance GTP-binding protein MxB) (Myxovirus resistance protein 2) (p78-related protein), MX2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	MX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.