



MYH11 Antibody, HRP conjugated

Product Code	CSB-PA015291LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35749
Immunogen	Recombinant Human Myosin-11 protein (1655-1797AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myosin-11 (Myosin heavy chain 11) (Myosin heavy chain, smooth muscle isoform) (SMMHC), MYH11, KIAA0866
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYH11