



MYH11 Antibody, Biotin conjugated

Product Code	CSB-PA015291LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35749
Immunogen	Recombinant Human Myosin-11 protein (1655-1797AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myosin-11 (Myosin heavy chain 11) (Myosin heavy chain, smooth muscle isoform) (SMMHC), MYH11, KIAA0866
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYH11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.