



# MYH11 Antibody, Biotin conjugated

|                     |  |
|---------------------|--|
| Product Code        | CSB-PA015291LD01HU   |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| Uniprot No.         | P35749   |
| Immunogen           | Recombinant Human Myosin-11 protein (1655-1797AA)  |
| Raised In           | Rabbit   |
| Species Reactivity  | Human  |
| Tested Applications | ELISA  |
| Form                | Liquid   |
| Conjugate           | Biotin   |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                       |
| Purification Method | >95%, Protein G purified   |
| Isotype             | IgG  |
| Clonality           | Polyclonal   |
| Alias               | Myosin-11 (Myosin heavy chain 11) (Myosin heavy chain, smooth muscle isoform) (SMMHC), MYH11, KIAA0866 |
| Immunogen Species   | Homo sapiens (Human)   |
| Research Area       | Signal Transduction  |
| Target Names        | MYH11  |