



NKX2-8 Antibody, HRP conjugated

Product Code	CSB-PA015847LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15522
Immunogen	Recombinant Human Homeobox protein Nkx-2.8 protein (1-160AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein Nkx-2.8 (Homeobox protein NK-2 homolog H), NKX2-8, NKX-2.8 NKX2G NKX2H
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NKX2-8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.