



NOMO1 Antibody, HRP conjugated

Product Code	CSB-PA015931LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15155
Immunogen	Recombinant Human Nodal modulator 1 protein (918-1044AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nodal modulator 1 (pM5 protein), NOMO1, PM5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NOMO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.