



NRCAM Antibody, Biotin conjugated

Product Code	CSB-PA016074LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92823
Immunogen	Recombinant Human Neuronal cell adhesion molecule protein (1194-1299AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuronal cell adhesion molecule (Nr-CAM) (Neuronal surface protein Bravo) (hBravo) (NgCAM-related cell adhesion molecule) (Ng-CAM-related), NRCAM, KIAA0343
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NRCAM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.