



PCCB Antibody, HRP conjugated

Product Code	CSB-PA017523LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05166
Immunogen	Recombinant Human Propionyl-CoA carboxylase beta chain, mitochondrial protein (200-314AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Propionyl-CoA carboxylase beta chain, mitochondrial (PCCase subunit beta) (EC 6.4.1.3) (Propanoyl-CoA:carbon dioxide ligase subunit beta), PCCB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PCCB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.