



# PCSK7 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA017646LB01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q16549
<b>Immunogen</b>	Recombinant Human Proprotein convertase subtilisin/kexin type 7 protein (189-362AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Proprotein convertase subtilisin/kexin type 7 (EC 3.4.21.-) (Lymphoma proprotein convertase) (Prohormone convertase 7) (Proprotein convertase 7) (PC7) (Proprotein convertase 8) (PC8) (hPC8) (Subtilisin/kexin-like protease PC7), PCSK7, LPC PC7 PC8 SPC7
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	PCSK7
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.